Image Tool - Bug #7023

Remember fit selection between acquisition types

09/17/2014 09:11 AM - John Diamond

Status:	Resolved	Start date:	09/17/2014
Priority:	Normal	Due date:	
Assignee:		% Done:	0%
Category:		Estimated time:	0.00 hour
Target version:		Spent time:	0.50 hour

Description

For a One-shot acquisition, it looks like you switch to looking at the image-tool fitting regardless of what's selected, which makes sense, but once we go back to Live Acquisition, it looks like the Image Tool fitting is still selected. I can get back to the front end fitting if I toggle the analysis radio buttons off front-end fitting and back again, but I was wondering if you could just detect which one is selected on live acquisition start.

History

#1 - 05/15/2015 02:09 PM - John Diamond

- Status changed from New to Resolved

Read the analysis source when a live acquisition is restarted and notify all interested listeners.

09/30/2020 1/1