

Main Injector Radiation Survey Application (PA4078, I128) - Support #7850

Delete all survey readings from 1/14/2015 survey that don't have fmr/h reading unit

02/12/2015 02:44 PM - Kyle Hazelwood

Status:	Closed	Start date:	02/12/2015
Priority:	Normal	Due date:	
Assignee:	Kyle Hazelwood	% Done:	0%
Category:		Estimated time:	0.00 hour
Target version:		Spent time:	0.50 hour
Description			
<p>Bruce tells me that survey readings should be done in "freeze" mode. When readings are properly done in freeze mode the unit is displayed as fmr/h. Improper readings have units of mr/h. The last survey done on 1/14/2015 has three readings with mr/h. I need to delete these entries.</p> <p>I should be able to edit the original CSV file and re-import the file. Matching timestamps override old data.</p>			

History

#1 - 02/24/2015 04:00 PM - Kyle Hazelwood

Done.

#2 - 02/24/2015 04:17 PM - Kyle Hazelwood

- Status changed from New to Closed